



DISPOSABLE EQUIPMENT COVERS

GC MEDICA ENTERPRISE LTD., (WUXI)
GC MEDICA ENTERPRISE LTD. (WUXI)

DISPOSABLE EQUIPMENT COVERS

GCMedica offers a range of sterile single use equipment covers to maintain instant sterile field for theatre equipment in operating room. Our products include camera cover, operating light handle cover, C-arm drape and camera drape. Easy to use design provides cost effective and time saving solutions.





Camera Cover

Designed with len and easy applying foam fixing ring. Available with square len and round len with different sizes for different camera.



Product code

Picture

Description

GC9006A01



Round len O.D 85mm, 10cm long

GC9006A02



Round len O.D 75mm, 8cm long

GC9006B01



Square len O.D 65mm, 10cm long



Light Handle Cover

Available with rigid light handle cover and flexible light handle cover.

Provide an instant sterile interface between non-sterile surgical lamp and the sterile surgeon and nurse.

We supply reusable metal adapters to standardize your lights to apply sterile light handle covers. Our adapters cover most operating room lights in market.









Product code	Picture	Description
GC9005B01		Flexible, green color with double rings
GC9005B02		Flexible, green color with rebar
GC9005B03		Flexible, green color
GC9005E01		Flexible, purple color
GC9005C01		Blue color
GC9005A01		Sponge fixing ring type, blue color
GC9005F01		Rigid, green color



Metal Adapter For Light Handle Cover

Standardize your lights to apply sterile light handle covers and quickly set up sterile field.
Various adapters fit most operating room lights in market.

Product code	Picture	Description	Applicable operating lights
GC1968A01		Fit for blue color light handle cover (GC9005C01)	MAQUET and STERIGIS AS lights
GC1968B01		Fit for blue color light handle cover (GC9005C01)	MAQUET and STERIGIS AS lights
GC1968C01		Fit for blue color light handle cover (GC9005C01)	MAQUET and STERIGIS AS lights
GC1968L01		Fit for green and purple flexible light handle cover (GC9005B01/ GC9005B02/9GC005E01/ GC9005S01)	XLED lights
GC1968L02		Fit for green and purple flexible light handle cover (GC9005B01/ GC9005B02/9GC005E01/ GC9005S01)	XLED lights
GC1968K01		Fit for green and purple flexible light handle cover (GC9005B01/ GC9005B02/9GC005E01/ GC9005S01)	XLED lights



Endoscopic Camera Drape

Endoscopic camera drape creates an instant sterile field around endoscopic camera heads and cabling. By eliminating the need for between-case sterilisation, these sleeve options can save money on costly repairs and prolong the life of equipment.





Quick-change Endoscopic Camera Drape

Buckle securing design allows quickly and easy scope change during procedures.

Closed buckle with high definition lens, acting as sterile barrier as well, offers camera full protection and extend its service life. 85mm broad ring reduces risk of contamination to the sterile field.



Product code

Picture

Description

GC1956A01



Size 150mm x 2.5m, Ring aperture 85mm



High Definition Endoscopic Camera Drape



High definition lens end ensure the HD image quality of camera. Closed end design offers camera full protection from contamination, acting as a reliable sterile barrier.

Allow quick and easy scope change during procedures without compromising sterility.

85mm broad ring reduces risk of contamination to the sterile field.

Product code

Picture

Description

GC1956B01




Size 160mm x 2.35m, Ring aperture 85mm



Ring Folding Endoscopic Camera Drapes

Ring style drapes ensure easy and efficient application without compromising sterility.

85mm broad ring reduces risk of contamination to the sterile field.




Product code

Picture

Description

GC1956C01



TPE elastomeric end with two tapes,
Drape film is contained within ring.

GC1956W01



Open tip with two tapes



C-Arm Drape

C-Arm drape is used for isolating X-Ray from the contamination source during the surgical operations. It quickly builds up a sterile field for patients and doctors during procedures. Clip fixing design for easy application.



Product code

GC9024M01

Description

Consist of 2 scope / fluoroscopy drapes
and one drape with 6 clips



DISPOSABLE EQUIPMENT COVERS

CONTACT US



LOUJIN INDUSTRIAL PARK, SHUOFANG,
WUXI, 214143, P.R.CHINA



0086-510-88880366



INFO@GCMEDICA.COM



WWW.GCMEDICA.COM
